



Dimer Dots UV Dose Verification: **OPERATING ROOM**

Date: _____ Facility: _____ Device: _____

Instructions for use: -Keep Dots away from ALL sources of light until ready to use -Place FLAT on Target Surface

#	Site	mJ/cm ²	✓	#	Site	mJ/cm ²	✓
1	Interior door handle main door to corridor			27	Mayo stand 2 tray top inclined edge #1		
2	Floor Corner 1			28	Mayo stand 2 tray top inclined edge #2		
3	Floor Corner 2			29	Mayo stand 2 tray top inclined edge #3		
4	Floor Corner 3			30	Mayo stand 2 tray top inclined edge #4		
5	Floor Corner 4			31	Mayo stand 2 tray undersurface horizontal center		
6	Floor 30 cm from mid patient table left			32	IV Pole elbow level circumferential #1		
7	Floor 30 cm from mid patient table right			33	IV Pole elbow level circumferential #2		
8	OR Patient table vertical edge left			34	IV Pole elbow level circumferential #3		
9	OR Patient table vertical edge right			35	IV Pole elbow level circumferential #4		
10	OR Patient table vertical edge at foot			36	Sitting stool top center horizontal		
11	OR Patient table vertical edge at head			37	Sitting stool vertical side 4 cm from top #1		
12	OR Patient table horizontal at head			38	Sitting stool vertical side 4 cm from top #2		
13	OR Patient table horizontal mid center			39	Sitting stool vertical side 4 cm from top #3		
14	OR Patient table horizontal at foot			40	Sitting stool vertical side 4 cm from top #4		
15	Back table vertical surface left front			41	Overhead light top surface 4 cm from edge		
16	Back table vertical surface right front			42	Overhead light bottom (lighted) surface 4 cm from edge		
17	Back table vertical surface left edge mid			43	Anesthesia keyboard mid		
18	Back table vertical surface right edge mid			44	Anesthesia machine vertical surface		
19	Back table top horizontal 20 cm from left edge mid			45	Anesthesia computer touch screen mid		
20	Back table top horizontal 20cm from right edge mid			46	Anesthesia chair mid horizontal surface		
21	Back table under surface horizontal 4 cm front edge			47	Bed adjustment control front surface		
22	Mayo stand 1 tray top inclined edge #1			48	Bed adjustment control rear surface		
23	Mayo stand 1 tray top inclined edge #2			49	Wall cabinet handle exterior door		
24	Mayo stand 1 tray top inclined edge #3			50	Wall cabinet interior door at elbow level away from hinge		
25	Mayo stand 1 tray top inclined edge #4						
26	Mayo stand 1 tray undersurface horizontal center						

TOTAL VERIFIED: _____ / 50





Dimer Dots UV Dose Verification: **OPERATING ROOM**

COMMENTS:

TOTAL VERIFIED: _____ / 50





Dimer Dots UV Dose Verification: **OPERATING ROOM**

Date: _____ Facility: _____ Device: _____

Instructions for use: -Keep Dots away from ALL sources of light until ready to use -Place FLAT on Target Surface

#	Site	mJ/cm ²	✓	#	Site	mJ/cm ²	✓
1				27			
2				28			
3				29			
4				30			
5				31			
6				32			
7				33			
8				34			
9				35			
10				36			
11				37			
12				38			
13				39			
14				40			
15				41			
16				42			
17				43			
18				44			
19				45			
20				46			
21				47			
22				48			
23				49			
24				50			
25							
26							

TOTAL VERIFIED: _____ / 50





Dimer Dots UV Dose Verification: **OPERATING ROOM**

TOTAL VERIFIED: _____ / **50**

