



## Dimer Dots UV Dose Verification: **PATIENT ROOM**

**Date:** \_\_\_\_\_ **Facility:** \_\_\_\_\_ **Device:** \_\_\_\_\_

*Instructions for use: -Keep Dots away from ALL sources of light until ready to use -Place FLAT on Target Surface*

#	Site	mJ/cm <sup>2</sup>	✓	#	Site	mJ/cm <sup>2</sup>	✓
1	Interior door handle			27	Visitor chair arm rest top horizontal near front edge		
2	Floor Corner 1			28	Assistance call button front		
3	Floor Corner 2			29	Assistance call button side		
4	Floor Corner 3			30	Vitals monitor controls or touch screen		
5	Floor Corner 4			31	TV or light control front surface		
6	Floor 30 cm from mid bed left			32	TV or light control rear surface		
7	Floor 30 cm from mid bed right			33	IV Pole elbow level circumferential #1		
8	Floor 30 cm from foot of bed			34	IV Pole elbow level circumferential #2		
9	Floor 30 cm from main door to corridor			35	IV Pole elbow level circumferential #3		
10	Bed rail mid toward patient left			36	IV Pole elbow level circumferential #4		
11	Bed rail mid toward patient right			37	Patient room light switch		
12	Bed footboard mid toward patient			38	Bathroom floor center		
13	Bed footboard top mid			39	Bathroom light switch		
14	Curtain elbow height leading edge toward bed			40	Bathroom door handle toward patient room		
15	Curtain elbow height leading edge away from bed			41	Bathroom door handle toward bathroom		
16	Curtain elbow height mid curtain toward bed			42	Bathroom shower curtain/door leading edge interior		
17	Curtain elbow height mid curtain away from bed			43	Bathroom shower curtain/door leading edge exterior		
18	Overbed table top center			44	Bathroom countertop near basin		
19	Overbed table vertical edge toward patient			45	Bathroom counter face near basin		
20	Overbed table interior floor of drawer			46	Bathroom basin bottom of basin near drain		
21	Overbed table vertical free edge			47	Bathroom basin valve handle top surface		
22	Overbed table undersurface near free edge			48	Bathroom toilet seat lower surface near front edge		
23	Night stand top center away from wall			49	Bathroom toilet seat upper surface near front edge		
24	Night stand front top drawer vertical face center			50	Bathroom toilet flush lever		
25	Night stand top drawer interior floor near drawer face						
26	Visitor chair mid seating surface horizontal						

TOTAL VERIFIED: \_\_\_\_\_ / 50





**COMMENTS:**





Dimer Dots UV Dose Verification: **PATIENT ROOM**

Date: \_\_\_\_\_ Facility: \_\_\_\_\_ Device: \_\_\_\_\_

Instructions for use: -Keep Dots away from ALL sources of light until ready to use -Place FLAT on Target Surface

#	Site	mJ/cm <sup>2</sup>	✓	#	Site	mJ/cm <sup>2</sup>	✓
1				27			
2				28			
3				29			
4				30			
5				31			
6				32			
7				33			
8				34			
9				35			
10				36			
11				37			
12				38			
13				39			
14				40			
15				41			
16				42			
17				43			
18				44			
19				45			
20				46			
21				47			
22				48			
23				49			
24				50			
25							
26							

TOTAL VERIFIED: \_\_\_\_\_ / 50





Dimer Dots UV Dose Verification: **PATIENT ROOM**

TOTAL VERIFIED: \_\_\_\_\_ / **50**

